

PFAC-MXK

Native PFAS Solution/Mixture

1,2mL

PFBA		Perfluoro-n-butanoic acid	2000ng/mL
PFPeA		Perfluoro-n-pentanoic acid	2000ng/mL
PFHxA		Perfluoro-n-hexanoic acid	2000ng/mL
PFHpA		Perfluoro-n-heptanoic acid	2000ng/mL
PFOA		Perfluoro-n-octanoic acid	2000ng/mL
PFNA		Perfluoro-n-nonanoic acid	2000ng/mL
HFPO-DA		2,3,3,3-Tetrafluoro-2-(1,1,2,2,3,3,3-heptafluoropropoxy)propanoic acid	2000ng/mL
L-PFBS		Potassium perfluoro-1-butanefulfonate	2000ng/mL
L-PFHxS		Sodium perfluoro-1-(1,2,3-13C3)hexanesulfonate	2000ng/mL
L-PFOS		Sodium perfluoro-1-octanesulfonate	2000ng/mL

NOTE: Listed concentrations for the perfluoroalkanesulfonates are reported as the salt.

MPFAC-MXK

Mass-Labelled PFAS Solution/Mixture

1,2mL

MPFBA		Perfluoro-n-(1,2,3,4-13C4)butanoic acid	2000ng/mL
M5PFPeA		Perfluoro-n-(1,2,3,4,5-13C5)pentanoic acid	2000ng/mL
M5PFHxA		Perfluoro-n-(1,2,3,4,6-13C5)hexanoic acid	2000ng/mL
M4PFHpA		Perfluoro-n-(1,2,3,4-13C4)heptanoic acid	2000ng/mL
M8PFOA		Perfluoro-n-(13C8)octanoic acid	2000ng/mL
M9PFNA		Perfluoro-n-(13C9)nonanoic acid	2000ng/mL
M3HFPO-DA		2,3,3,3-Tetrafluoro-2-(1,1,2,2,3,3,3-heptafluoropropoxy)(13C3)propanoic acid	2000ng/mL
M3PFBS		Sodium perfluoro-1-(2,3,4-13C3)butanesulfonate	2000ng/mL
M3PFHxS		Sodium perfluoro-1-(1,2,3-13C3)hexanesulfonate	2000ng/mL
M8PFOS		Sodium perfluoro-(13C8)octanesulfonate	2000ng/mL

NOTE: Listed concentrations for the perfluoroalkanesulfonates are reported as the salt.