




## **NEW PRODUCT**


### **7H-PFHpA**

In 2025, the European Committee for Standardization (CEN) approved EN 17681-1:2025 for the analysis of PFAS in textiles and textile products. Included in the updated list of target compounds is 7H-perfluoroheptanoic acid (**7H-PFHpA** or HPFHpA). Additionally, **7H-PFHpA** is listed as an analyte in the ISO/DIS 25652:2025 method for the analysis of PFAS in sediment, soil, sludge and waste. In response to these new methods, **Wellington** is pleased to announce the release of a certified reference standard for **7H-PFHpA**.



**7H-PFHpA**

Product Code	Description (methanol solution)	Units	Conc
<b>NEW</b> 7H-PFHpA	7H-Perfluoroheptanoic acid	1.2 mL	50.0 µg/mL



Full scan chromatogram, extracted mass spectrum, and selected MRM transitions of **7H-PFHpA**.

*Please contact your local distributor or [info@well-labs.com](mailto:info@well-labs.com) for pricing and delivery.*

*Visit our website ([www.well-labs.com](http://www.well-labs.com)) for a complete listing of new products.*