

有機フッ素化合物 標準液
Perfluorinated Compounds(PFCs)

2021



 WELLINGTON
LABORATORIES

Wellington Laboratories Japan Inc.
www.well-labs.co.jp

PFC-CVS-C

Catalogue Number	Product (methanol solution)	Qty/Conc				
用 PFC-CVS-C	PFC-C-CS5 Calibration Solution CS1 - CS5	1 kit (5 ampoules)				
PFC-C-CS1	CS1	200µl				
PFC-C-CS2	CS2	200µl				
PFC-C-CS3	CS3	200µl				
PFC-C-CS4	CS4	200µl				
PFC-C-CS5	CS5	200µl				
Native PFCs		PFC-C-CS1 (ng/ml)	PFC-C-CS2 (ng/ml)	PFC-C-CS3 (ng/ml)	PFC-C-CS4 (ng/ml)	PFC-C-CS5 (ng/ml)
Perfluoro-n-butanoic acid	PFBA	2.0	10	50	200	1000
Perfluoro-n-pentanoic acid	PPPeA	2.0	10	50	200	1000
Perfluoro-n-hexanoic acid	PFHxA	2.0	10	50	200	1000
Perfluoro-n-heptanoic acid	PFHpA	2.0	10	50	200	1000
Perfluoro-n-octanoic acid	PFOA	2.0	10	50	200	1000
Perfluoro-n-nonanoic acid	PFNA	2.0	10	50	200	1000
Perfluoro-n-decanoic acid	PFDA	2.0	10	50	200	1000
Perfluoro-n-undecanoic acid	PFUdA	2.0	10	50	200	1000
Perfluoro-n-dodecanoic acid	PFDoA	2.0	10	50	200	1000
Perfluoro-n-tridecanoic acid	PTFDA	2.0	10	50	200	1000
Perfluoro-n-tetradecanoic acid	PFTeDA	2.0	10	50	200	1000
Perfluoro-n-hexadecanoic acid	PFHxDA	2.0	10	50	200	1000
Perfluoro-n-octadecanoic acid	PFODA	2.0	10	50	200	1000
Potassium perfluoro-1-butanesulfonate	L-PFBS	2.0	10	50	200	1000
Sodium perfluoro-1-pentanesulfonate	L-PFPeS	2.0	10	50	200	1000
Sodium perfluoro-1-hexanesulfonate	L-PFHxS	2.0	10	50	200	1000
Sodium perfluoro-1-heptanesulfonate	L-PFHxS	2.0	10	50	200	1000
Sodium perfluoro-1-octanesulfonate	L-PFOS	2.0	10	50	200	1000
Sodium perfluoro-1-nonanesulfonate	L-PFNS	2.0	10	50	200	1000
Sodium perfluoro-1-decanesulfonate	L-PFDS	2.0	10	50	200	1000
Sodium perfluoro-1-dodecanesulfonate	L-PFDs	2.0	10	50	200	1000
Mass-Labelled PFCs Extraction Standards						
Perfluoro-n-[1,2,3,4- ¹³ C ₄]butanoic acid	MPFBA	50	50	50	50	50
Perfluoro-n-[1,2,3,4,5- ¹³ C ₅]pentanoic acid	M5PPeA	50	50	50	50	50
Perfluoro-n-[1,2,3,4,6- ¹³ C ₅]hexanoic acid	M5PFHxA	50	50	50	50	50
Perfluoro-n-[1,2,3,4- ¹³ C ₄]heptanoic acid	M4PFHpA	50	50	50	50	50
Perfluoro-n-[¹³ C ₈]octanoic acid	M8PFOA	50	50	50	50	50
Perfluoro-n-[¹³ C ₉]nonanoic acid	M9PFNA	50	50	50	50	50
Perfluoro-n-[1,2,3,4,5,6- ¹³ C ₆]decanoic acid	M6PFDA	50	50	50	50	50
Perfluoro-n-[1,2,3,4,5,6,7- ¹³ C ₇]undecanoic acid	M7PFUdA	50	50	50	50	50
Perfluoro-n-[1,2- ¹³ C ₂]dodecanoic acid	MPFDoA	50	50	50	50	50
Perfluoro-n-[1,2- ¹³ C ₂]tetradecanoic acid	M2PFTeDA	50	50	50	50	50
Sodium perfluoro-1-[2,3,4- ¹³ C ₃]butanesulfonate	M3PFBS	50	50	50	50	50
Sodium perfluoro-1-[1,2,3- ¹³ C ₃]hexanesulfonate	M3PFHxS	50	50	50	50	50
Sodium perfluoro-[¹³ C ₈]octanesulfonate	M8PFOS	50	50	50	50	50
Mass-Labelled PFCs Injection Standards						
Perfluoro-n-[1,2,3,4- ¹³ C ₃]butanoic acid	M3PFBA	50	50	50	50	50
Perfluoro-n-[1,2- ¹³ C ₂]octanoic acid	M2PFOA	50	50	50	50	50
Perfluoro-n-[1,2- ¹³ C ₂]decanoic acid	MPFDA	50	50	50	50	50
Sodium perfluoro-1-[1,2,3,4- ¹³ C ₄]octanesulfonate	MPFOS	50	50	50	50	50

NOTE: Listed concentrations are reported as the salt.

用 第一種特定化學物質

PFC-CVS-C

Catalogue Number	Product (methanol solution)	Qty/Conc		
用 MPFAC-C-ES	Mass-Labelled PFC Extraction Standards Soluiton	—	1.2ml	
用 MPFAC-C-IS	Mass-Labelled PFC Injection Standards Soluiton	—	1.2ml	
用 PFAC-MXC	Native PFC Stock Soluiton	—	1.2ml	
Native PFCs		MPFAC-C-ES (ng/ml)	MPFAC-C-IS (ng/ml)	PFAC-MXC (ng/ml)
Perfluoro-n-butanoic acid	PFBA	—	—	2000
Perfluoro-n-pentanoic acid	PPPeA	—	—	2000
Perfluoro-n-hexanoic acid	PFHxA	—	—	2000
Perfluoro-n-heptanoic acid	PFHpA	—	—	2000
Perfluoro-n-octanoic acid	PFOA	—	—	2000
Perfluoro-n-nonanoic acid	PFNA	—	—	2000
Perfluoro-n-decanoic acid	PFDA	—	—	2000
Perfluoro-n-undecanoic acid	PFUdA	—	—	2000
Perfluoro-n-dodecanoic acid	PFDoA	—	—	2000
Perfluoro-n-tridecanoic acid	PFTrDA	—	—	2000
Perfluoro-n-tetradecanoic acid	PFTeDA	—	—	2000
Perfluoro-n-hexadecanoic acid	PFHxDA	—	—	2000
Perfluoro-n-octadecanoic acid	PFODA	—	—	2000
Potassium perfluoro-1-butanesulfonate	L-PFBS	—	—	2000
Sodium perfluoro-1-pentanesulfonate	L-PFPeS	—	—	2000
Sodium perfluoro-1-hexanesulfonate	L-PFHxS	—	—	2000
Sodium perfluoro-1-heptanesulfonate	L-PFHps	—	—	2000
Sodium perfluoro-1-octanesulfonate	L-PFOS	—	—	2000
Sodium perfluoro-1-nonanesulfonate	L-PFNS	—	—	2000
Sodium perfluoro-1-decanesulfonate	L-PFDS	—	—	2000
Sodium perfluoro-1-dodecanesulfonate	L-PFDsO	—	—	2000
Mass-Labelled PFCs Extraction Standards				
Perfluoro-n-[1,2,3,4- ¹³ C ₄]butanoic acid	MPFBA	2000	—	—
Perfluoro-n-[1,2,3,4,5- ¹³ C ₅]pentanoic acid	M5PFPeA	2000	—	—
Perfluoro-n-[1,2,3,4,6- ¹³ C ₅]hexanoic acid	M5PFHxA	2000	—	—
Perfluoro-n-[1,2,3,4- ¹³ C ₄]heptanoic acid	M4PFHpA	2000	—	—
Perfluoro-n-[¹³ C ₈]octanoic acid	M8PFOA	2000	—	—
Perfluoro-n-[¹³ C ₉]nonanoic acid	M9PFNA	2000	—	—
Perfluoro-n-[1,2,3,4,5,6- ¹³ C ₆]decanoic acid	M6PFDA	2000	—	—
Perfluoro-n-[1,2,3,4,5,6,7- ¹³ C ₇]undecanoic acid	M7PFUdA	2000	—	—
Perfluoro-n-[1,2- ¹³ C ₂]dodecanoic acid	MPFDoA	2000	—	—
Perfluoro-n-[1,2- ¹³ C ₂]tetradecanoic acid	M2PFTeDA	2000	—	—
Sodium perfluoro-1-[2,3,4- ¹³ C ₃]butanesulfonate	M3PFBS	2000	—	—
Sodium perfluoro-1-[1,2,3- ¹³ C ₃]hexanesulfonate	M3PFHxS	2000	—	—
Sodium perfluoro-[¹³ C ₈]octanesulfonate	M8PFOS	2000	—	—
Mass-Labelled PFCs Injection Standards				
Perfluoro-n-[1,2,3,4- ¹³ C ₃]butanoic acid	M3PFBA	—	2000	—
Perfluoro-n-[1,2- ¹³ C ₂]octanoic acid	M2PFOA	—	2000	—
Perfluoro-n-[1,2- ¹³ C ₂]decanoic acid	MPFDA	—	2000	—
Sodium perfluoro-1-[1,2,3,4- ¹³ C ₄]octanesulfonate	MPFOS	—	2000	—

NOTE: Listed concentrations are reported as the salt.

NATIVE LINEAR PERFLUOROALKYLSULFONATES (PFSAs)

Catalogue Number	Product (methanol solution)	Qty/Conc	
L-PFPPrS	Sodium perfluoro-1-propanesulfonate	1.2 ml	50 µg/ml
L-PFBS	Potassium perfluoro-1-butanesulfonate	1.2 ml	50 µg/ml
L-PFPeS	Sodium perfluoro-1-pentanesulfonate	1.2 ml	50 µg/ml
L-PFHxS	Sodium perfluoro-1-hexanesulfonate	1.2 ml	50 µg/ml
L-PFHpS	Sodium perfluoro-1-heptanesulfonate	1.2 ml	50 µg/ml
用 L-PFOS	Sodium perfluoro-1-octanesulfonate	1.2 ml	50 µg/ml
用 L-PFOSK	Potassium perfluoro-1-octanesulfonate	1.2 ml	50 µg/ml
L-PFNS	Sodium perfluoro-1-nonanesulfonate	1.2 ml	50 µg/ml
L-PFDS	Sodium perfluoro-1-decanesulfonate	1.2 ml	50 µg/ml
NEW L-PFUDs	Sodium perfluoro-1-undecanesulfonate	1.2 ml	50 µg/ml
L-PFDoS	Sodium perfluoro-1-dodecanesulfonate	1.2 ml	50 µg/ml
NEW L-PFTTrDS	Sodium perfluoro-1-tridecanesulfonate	1.2 ml	50 µg/ml

NOTE: Listed concentrations are reported as the salt.

NATIVE PERFLUOROALKYLSULFONATES: SOLUTION/MIXTURE

Catalogue Number	Product (methanol solution)	Qty/Conc	
用 PFS-MXA	Native PFAS Solution/Mixture	1.2 ml	
Potassium perfluoro-1-butanesulfonate		2 µg/ml	
Sodium perfluoro-1-hexanesulfonate		2 µg/ml	
Sodium perfluoro-1-heptanesulfonate		2 µg/ml	
Sodium perfluoro-1-octanesulfonate		2 µg/ml	
Sodium perfluoro-1-decanesulfonate		2 µg/ml	

NOTE: Listed concentrations are reported as the salt.

NATIVE CHLORO-PERFLUOROALKYLSULFONATE

Catalogue Number	Product (methanol solution)	Qty/Conc	
8CI-PFOS	Sodium 8-chloroperfluoro-1-octanesulfonate	1.2 ml	50 µg/ml

NOTE: Listed concentrations are reported as the salt.

NATIVE PERFLUOROETHYLCYCLOHEXANESULFONATE

Catalogue Number	Product (methanol solution)	Qty/Conc	
PFECHS	Perfluoro-4-ethylcyclohexanesulfonate	1.2 ml	50 µg/ml






NOTE: Listed concentrations are reported as the salt.

NATIVE BRANCHED PERFLUOROALKYLSULFONATES

Catalogue Number	Product (methanol solution)	Qty/Conc	
br-PFHxSK	L-PFHxS with branched isomers (Potassium Salt)	1.2 ml	50 µg/ml
用 br-PFOSK	L-PFOSK with branched isomers	1.2 ml	50 µg/ml
用 T-PFOS	Potassium perfluorooctanesulfonate (Technical Grade)	1.2 ml	50 µg/ml
NaP3MHpS	Sodium perfluoro-3-methylheptanesulfonate	1.2 ml	50 µg/ml
NaP6MHpS	Sodium perfluoro-6-methylheptanesulfonate	1.2 ml	50 µg/ml
ipPFNS	Sodium perfluoro-7-methyloctanesulfonate	1.2 ml	50 µg/ml

NOTE: Listed concentrations are reported as the salt.

MASS-LABELLED PERFLUOROALKYLSULFONATES

Catalogue Number	Product
M3PFBS	 <p>Sodium perfluoro-1-[2,3,4-¹³C₃]butanesulfonate 1.2 ml; 50 µg/ml (± 2.5 µg/ml); in methanol >99% linear; >99% ¹³C₃</p>
MPFHxS	 <p>Sodium perfluoro-1-hexane[¹⁸O₂]sulfonate 1.2 ml; 50 µg/ml (± 2.5 µg/ml); in methanol >99% linear; >94% ¹⁸O₂</p>
M3PFHxS	 <p>Sodium perfluoro-1-[1,2,3-¹³C₃]hexanesulfonate 1.2 ml; 50 µg/ml (± 2.5 µg/ml); in methanol >99% linear; >99% ¹³C₃</p>
用 MPFOS	 <p>Sodium perfluoro-1-[1,2,3,4-¹³C₄]- octanesulfonate 1.2 ml; 50 µg/ml (± 2.5 µg/ml); in methanol >99% linear; >99% ¹³C₄</p>
用 M8PFOS	 <p>Sodium perfluoro-[¹³C₈]octanesulfonate 1.2 ml; 50 µg/ml (± 2.5 µg/ml); in methanol >99% linear; >99% ¹³C₈</p>

NOTE: Listed concentrations are reported as the salt.

NATIVE LINEAR PERFLUOROALKYLCARBOXYLIC ACIDS (PFCAs)

Catalogue Number	Product (methanol solution)	Qty/Conc	
PFBA	Perfluoro-n-butanoic acid	1.2 ml	50 µg/ml
PPPeA	Perfluoro-n-pentanoic acid	1.2 ml	50 µg/ml
PFHxA	Perfluoro-n-hexanoic acid	1.2 ml	50 µg/ml
PFHpA	Perfluoro-n-heptanoic acid	1.2 ml	50 µg/ml
PFOA	Perfluoro-n-octanoic acid	1.2 ml	50 µg/ml
PFNA	Perfluoro-n-nonanoic acid	1.2 ml	50 µg/ml
PFDA	Perfluoro-n-decanoic acid	1.2 ml	50 µg/ml
PFUdA	Perfluoro-n-undecanoic acid	1.2 ml	50 µg/ml
PFDoA	Perfluoro-n-dodecanoic acid	1.2 ml	50 µg/ml
PFTrDA	Perfluoro-n-tridecanoic acid	1.2 ml	50 µg/ml
PFTeDA	Perfluoro-n-tetradecanoic acid	1.2 ml	50 µg/ml
PFHxDA	Perfluoro-n-hexadecanoic acid	1.2 ml	50 µg/ml
PFODA	Perfluoro-n-octadecanoic acid	1.2 ml	50 µg/ml

 第一種特定化學物質






NATIVE PERFLUOROALKYLCARBOXYLIC ACIDS: SOLUTION/MIXTURE

Catalogue Number	Product (methanol solution)	Qty/Conc
PFC-MXA	Native PFCA Solution/Mixture	1.2 ml
Perfluoro-n-butanoic acid		2 µg/ml
Perfluoro-n-pentanoic acid		2 µg/ml
Perfluoro-n-hexanoic acid		2 µg/ml
Perfluoro-n-heptanoic acid		2 µg/ml
Perfluoro-n-octanoic acid		2 µg/ml
Perfluoro-n-nonanoic acid		2 µg/ml
Perfluoro-n-decanoic acid		2 µg/ml
Perfluoro-n-undecanoic acid		2 µg/ml
Perfluoro-n-dodecanoic acid		2 µg/ml
Perfluoro-n-tridecanoic acid		2 µg/ml
Perfluoro-n-tetradecanoic acid		2 µg/ml

NATIVE BRANCHED PERFLUOROALKYLCARBOXYLIC ACIDS

Catalogue Number	Product (methanol solution)	Qty/Conc
T-PFOA	Technical Ammonium Perfluorooctanoate	1.2 ml 50 µg/ml
P3MHPa	Perfluoro-3-methylheptanoic acid	1.2 ml 50 µg/ml
P4MOA	Perfluoro-4-methyloctanoic acid	1.2 ml 50 µg/ml
ipPFNA	Perfluoro-7-methyloctanoic acid	1.2 ml 50 µg/ml
P355TMHxA	Perfluoro-3,5,5-trimethylhexanoic acid	1.2 ml 50 µg/ml
P37DMOA	Perfluoro-3,7-dimethyloctanoic acid	1.2 ml 50 µg/ml

MASS-LABELLED PERFLUOROALKYLCARBOXYLIC ACIDS

Catalogue Number	Product
M3PFBA	 Perfluoro-n-[2,3,4- ¹³ C ₃]butanoic acid 1.2 ml; 50 µg/ml ($\pm 2.5 \mu\text{g/ml}$); in methanol >99% linear; >99% 2,3,4- ¹³ C ₃
MPFBA	 Perfluoro-n-[1,2,3,4- ¹³ C ₄]butanoic acid 1.2 ml; 50 µg/ml ($\pm 2.5 \mu\text{g/ml}$); in methanol >99% linear; >99% ¹³ C ₄
M3PFPeA	 Perfluoro-n-[3,4,5- ¹³ C ₃]pentanoic acid 1.2 ml; 50 µg/ml ($\pm 2.5 \mu\text{g/ml}$); in methanol >99% linear; >99% 3,4,5- ¹³ C ₃
M5PFPeA	 Perfluoro-n-[1,2,3,4,5- ¹³ C ₅]pentanoic acid 1.2 ml; 50 µg/ml ($\pm 2.5 \mu\text{g/ml}$); in methanol >99% linear; >99% ¹³ C ₅
MPFHxA	 Perfluoro-n-[1,2- ¹³ C ₂]hexanoic acid 1.2 ml; 50 µg/ml ($\pm 2.5 \mu\text{g/ml}$); in methanol >99% linear; >99% ¹³ C ₂

MASS-LABELLED PERFLUOROALKYLCARBOXYLIC ACIDS

Catalogue Number	Product
M5PFHxA	Perfluoro-n-[1,2,3,4,6- ¹³ C ₅]hexanoic acid 1.2 ml; 50 µg/ml (\pm 2.5 µg/ml); in methanol >99% linear; >99% ¹³ C ₅
M4PFHpA	Perfluoro-n-[1,2,3,4,- ¹³ C ₄]heptanoic acid 1.2 ml; 50 µg/ml (\pm 2.5 µg/ml); in methanol >99% linear; >99% ¹³ C ₄
M2PFOA	Perfluoro-n-[1,2- ¹³ C ₂]octanoic acid 1.2 ml; 50 µg/ml (\pm 2.5 µg/ml); in methanol >99% linear; >99% ¹³ C ₂
MPFOA	Perfluoro-n-[1,2,3,4- ¹³ C ₄]octanoic acid 1.2 ml; 50 µg/ml (\pm 2.5 µg/ml); in methanol >99% linear; >99% ¹³ C ₄
M8PFOA	Perfluoro-n-[¹³ C ₈]octanoic acid 1.2 ml; 50 µg/ml (\pm 2.5 µg/ml); in methanol >99% linear; 97.9% ¹³ C ₈ and 2.1% ¹³ C ₄
MPFNA	Perfluoro-n-[1,2,3,4,5- ¹³ C ₅]nonanoic acid 1.2 ml; 50 µg/ml (\pm 2.5 µg/ml); in methanol >99% linear; >99% ¹³ C ₅
M9PFNA	Perfluoro-n-[¹³ C ₉]nonanoic acid 1.2 ml; 50 µg/ml (\pm 2.5 µg/ml); in methanol >99% linear; >99% ¹³ C ₉
MPFDA	Perfluoro-n-[1,2- ¹³ C ₂]decanoic acid 1.2 ml; 50 µg/ml (\pm 2.5 µg/ml); in methanol >99% linear; >99% ¹³ C ₂
MPFDA-A	Perfluoro-n-[1,2- ¹³ C ₂]decanoic acid 1.2 ml; 50 µg/ml (\pm 2.5 µg/ml); in acetonitrile >99% linear; >99% ¹³ C ₂
M6PFDA	Perfluoro-n-[1,2,3,4,5,6- ¹³ C ₆]decanoic acid 1.2 ml; 50 µg/ml (\pm 2.5 µg/ml); in methanol >99% linear; >99% ¹³ C ₆
MPFUdA	Perfluoro-n-[1,2- ¹³ C ₂]undecanoic acid 1.2 ml; 50 µg/ml (\pm 2.5 µg/ml); in methanol >99% linear; >99% ¹³ C ₂
M7PFUdA	Perfluoro-n-[1,2,3,4,5,6,7- ¹³ C ₇]undecanoic acid 1.2 ml; 50 µg/ml (\pm 2.5 µg/ml); in methanol >99% linear; >99% ¹³ C ₇
MPFDa	Perfluoro-n-[1,2- ¹³ C ₂]dodecanoic acid 1.2 ml; 50 µg/ml (\pm 2.5 µg/ml); in methanol >99% linear; >99% ¹³ C ₂
M2PFTeDA	Perfluoro-n-[1,2- ¹³ C ₂]tetradecanoic acid 1.2 ml; 50 µg/ml (\pm 2.5 µg/ml); in methanol >99% linear; >99% ¹³ C ₂
M2PFHxDA	Perfluoro-n-[1,2- ¹³ C ₂]hexadecanoic acid 1.2 ml; 50 µg/ml (\pm 2.5 µg/ml); in methanol >99% linear; >99% ¹³ C ₂

MIXED NATIVE PFCAs AND PFASs: SOLUTION/MIXTURE

Catalogue Number	Product (methanol solution)	Qty/Conc
用 PFAC-MXA	Native PFCAs and PFASs Solution/Mixture	1.2 ml
Perfluoro-n-butanoic acid	PFBA	5 µg/ml
Perfluoro-n-pentanoic acid	PPPeA	5 µg/ml
Perfluoro-n-hexanoic acid	PFHxA	5 µg/ml
Perfluoro-n-heptanoic acid	PFHpA	5 µg/ml
Perfluoro-n-octanoic acid	PFOA	5 µg/ml
Perfluoro-n-nonanoic acid	PFNA	5 µg/ml
Perfluoro-n-decanoic acid	PFDA	5 µg/ml
Potassium perfluoro-1-butanesulfonate	L-PFBS	5 µg/ml
Sodium perfluoro-1-hexanesulfonate	L-PFHxS	5 µg/ml
Sodium perfluoro-1-octanesulfonate	L-PFOS	5 µg/ml
用 PFAC-MXB	Native PFCAs and PFASs Solution/Mixture	1.2 ml
Perfluoro-n-butanoic acid	PFBA	2 µg/ml
Perfluoro-n-pentanoic acid	PPPeA	2 µg/ml
Perfluoro-n-hexanoic acid	PFHxA	2 µg/ml
Perfluoro-n-heptanoic acid	PFHpA	2 µg/ml
Perfluoro-n-octanoic acid	PFOA	2 µg/ml
Perfluoro-n-nonanoic acid	PFNA	2 µg/ml
Perfluoro-n-decanoic acid	PFDA	2 µg/ml
Perfluoro-n-undecanoic acid	PFUdA	2 µg/ml
Perfluoro-n-dodecanoic acid	PFDoA	2 µg/ml
Perfluoro-n-tridecanoic acid	PFTrDA	2 µg/ml
Perfluoro-n-tetradecanoic acid	PFTeDA	2 µg/ml
Perfluoro-n-hexadecanoic acid	PFHxDA	2 µg/ml
Perfluoro-n-octadecanoic acid	PFODA	2 µg/ml
Potassium perfluoro-1-butanesulfonate	L-PFBS	2 µg/ml
Sodium perfluoro-1-hexanesulfonate	L-PFHxA	2 µg/ml
Sodium perfluoro-1-octanesulfonate	L-PFOS	2 µg/ml
Sodium perfluoro-1-decanesulfonate	L-PFDS	2 µg/ml
NEW 用 EU-5813-NSS	5813/20 PFAS Native Solution/Mixture	1.2 ml
Perfluoro-n-butanoic acid	PFBA	2 µg/ml
Perfluoro-n-pentanoic acid	PPPeA	2 µg/ml
Perfluoro-n-hexanoic acid	PFHxA	2 µg/ml
Perfluoro-n-heptanoic acid	PFHxA	2 µg/ml
Perfluoro-n-octanoic acid	PFHpA	2 µg/ml
Perfluoro-n-nonanoic acid	PFOA	2 µg/ml
Perfluoro-n-decanoic acid	PFNA	2 µg/ml
Perfluoro-n-undecanoic acid	PFDA	2 µg/ml
Perfluoro-n-dodecanoic acid	PFDoA	2 µg/ml
Perfluoro-n-tridecanoic acid	PFTrDA	2 µg/ml
Potassium perfluoro-1-butanesulfonate	L-PFBS	2 µg/ml
Sodium perfluoro-1-pentanesulfonate	L-PFPeS	2 µg/ml
Sodium perfluoro-1-hexanesulfonate	L-PFHxS	2 µg/ml
Sodium perfluoro-1-heptanesulfonate	L-PFHpS	2 µg/ml
Sodium perfluoro-1-octanesulfonate	L-PFOS	2 µg/ml
Sodium perfluoro-1-nonanesulfonate	L-PFNS	2 µg/ml
Sodium perfluoro-1-decanesulfonate	L-PFDS	2 µg/ml
Sodium perfluoro-1-undecanesulfonate	L-PFUdS	2 µg/ml
Sodium perfluoro-1-dodecanesulfonate	L-PFDoS	2 µg/ml
Sodium perfluoro-1-tridecanesulfonate	L-PFTrDS	2 µg/ml

MIXED MASS-LABELLED PFCAs AND PFASs: SOLUTION/MIXTURE

Catalogue Number	Product (methanol solution)	Qty/Conc
用 MPFAC-MXA	Mass-Labelled PFCAs and PFASs Solution/Mixture	1.2 ml
Perfluoro-n-[1,2,3,4- ¹³ C ₄]butanoic acid	MPFBA	2 µg/ml
Perfluoro-n-[1,2- ¹³ C ₂]hexanoic acid	MPFHxA	2 µg/ml
Perfluoro-n-[1,2,3,4- ¹³ C ₄]octanoic acid	MPFOA	2 µg/ml
Perfluoro-n-[1,2,3,4,5- ¹³ C ₅]nonanoic acid	MPFNA	2 µg/ml
Perfluoro-n-[1,2- ¹³ C ₂]decanoic acid	MPFDA	2 µg/ml
Perfluoro-n-[1,2- ¹³ C ₂]undecanoic acid	MPFUdA	2 µg/ml
Perfluoro-n-[1,2- ¹³ C ₂]dodecanoic acid	MPFDoA	2 µg/ml
Sodium perfluoro-1-hexane[¹⁸ O ₂]sulfonate	MPFHxS	2 µg/ml
Sodium perfluoro-1-[1,2,3,4- ¹³ C ₄]octanesulfonate	MPFOS	2 µg/ml

NOTE: Listed concentrations for the perfluoroalkylsulfonates are reported as the salt.

MIXED NATIVE PFAS: SOLUTION/MIXTURES

Catalogue Number	Product (methanol/4% isopropanol,/ <1% water)	Qty/Conc
用 PFAC-24PAR	Native PFAS Precision and Recovery Standard Solution	1.2 ml
Perfluoro-n-butanoic acid	PFBA	2000
Perfluoro-n-pentanoic acid	PPPeA	2000
Perfluoro-n-hexanoic acid	PFHxA	2000
Perfluoro-n-heptanoic acid	PFHpA	2000
Perfluoro-n-octanoic acid	PFOA	2000
Perfluoro-n-nonanoic acid	PFNA	2000
Perfluoro-n-decanoic acid	PFDA	2000
Perfluoro-n-undecanoic acid	PFUdA	2000
Perfluoro-n-dodecanoic acid	PFDoA	2000
Perfluoro-n-tridecanoic acid	PFTrDA	2000
Perfluoro-n-tetradecanoic acid	PFTeDA	2000
Perfluoro-1-octanesulfonamide	FOSA	2000
N-methylperfluoro-1-octanesulfonamidoacetic acid	N-MeFOSAA	2000
N-ethylperfluoro-1-octanesulfonamidoacetic acid	N-EtFOSAA	2000
Potassium perfluoro-1-butanesulfonate	L-PFBS	2000
Sodium perfluoro-1-pentanesulfonate	L-PFPeS	2000
Potassium perfluorohexanenesulfonate	PFHxSK: linear isomer	1620
	PFHxSK: Σ branched isomers	378
Sodium perfluoro-1-heptanesulfonate	L-PFHpS	2000
Potassium perfluoroctanenesulfonate	PFOSK: linear isomer	1580
	PFOSK: Σ branched isomers	422
Sodium perfluoro-1-nananesulfonate	L-PFNS	2000
Sodium perfluoro-1-decanesulfonate	L-PFDS	2000
Sodium 1H,1H,2H,2H-perfluorohexane sulfonate	4:2FTS	2000
Sodium 1H,1H,2H,2H-perfluorooctane sulfonate	6:2FTS	2000
Sodium 1H,1H,2H,2H-perfluorodecane sulfonate	8:2FTS	2000

NOTE: Listed concentrations are reported as the salt.

用 第一種特定化學物質

MIXED NATIVE PFAS: SOLUTION/MIXTURES

Catalogue Number	Product (methanol/6% isopropanol,/ <1% water)	Qty/Conc
用 PFAC-30PAR	Native PFAS Precision and Recovery Standard Solution	1.2 ml
Perfluoro-n-butanoic acid	PFBA	1000
Perfluoro-n-pentanoic acid	PPeA	1000
Perfluoro-n-hexanoic acid	PFHxA	1000
Perfluoro-n-heptanoic acid	PFHpA	1000
Perfluoro-n-octanoic acid	PFOA	1000
Perfluoro-n-nonanoic acid	PFNA	1000
Perfluoro-n-decanoic acid	PFDA	1000
Perfluoro-n-undecanoic acid	PFUdA	1000
Perfluoro-n-dodecanoic acid	PFDoA	1000
Perfluoro-n-tridecanoic acid	PFTrDA	1000
Perfluoro-n-tetradecanoic acid	PFTeDA	1000
Perfluoro-1-butanesulfonamide	FBSA	1000
Perfluoro-1-hexanesulfonamide	FHxSA	1000
Perfluoro-1-octanesulfonamide	FOSA	1000
2,3,3,3-Tetrafluoro-2-(1,1,2,2,3,3,3-heptafluoropropoxy)propanoic acid	HFPO-DA	1000
N-methylperfluoro-1-octanesulfonamidoacetic acid	N-MeFOSAA: linear isomer	760
	N-MeFOSAA: Σ branched isomers	240
N-ethylperfluoro-1-octanesulfonamidoacetic acid	N-EtFOSAA: linear isomer	775
	N-EtFOSAA: Σ branched isomers	225
Potassium perfluoro-1-butanesulfonate	L-PFBS	1000
Sodium perfluoro-1-pentanesulfonate	L-PPeS	1000
Potassium perfluorohexanenesulfonate	PFHxSK: linear isomer	811
	PFHxSK: Σ branched isomers	189
Sodium perfluoro-1-heptanesulfonate	L-PFHxS	1000
Potassium perfluoroctanenesulfonate	PFOSK: linear isomer	788
	PFOSK: Σ branched isomers	211
Sodium perfluoro-1-nonanesulfonate	L-PFNS	1000
Sodium perfluoro-1-decanesulfonate	L-PFDS	1000
Sodium 1H,1H,2H,2H-perfluorohexane sulfonate	4:2FTS	1000
Sodium 1H,1H,2H,2H-perfluorooctane sulfonate	6:2FTS	1000
Sodium 1H,1H,2H,2H-perfluorodecane sulfonate	8:2FTS	1000
Sodium dodecafluoro-3H-4,8-dioxanonanoate	NaDONA	1000
Potassium 9-chlorohexadecafluoro-3-oxanonane-1-sulfonate	9Cl-PF3ONS	1000
Potassium 11-chloroeicosafafluoro-3-oxaundecane-1-sulfonate	11Cl-PF3OUdS	1000

NOTE: Listed concentrations are reported as the salt.

用 第一種特定化学物質

MIXED NATIVE PFAS: SOLUTION/MIXTURES

Catalogue Number	Product (methanol <1% water)	Qty/Conc
用 EPA-537PDS-R1	EPA Method 537 Analyte Primary Dilution Standard (branched/linear mix)	1.2 ml
Perfluoro-n-hexanoic acid	PFHxA	2000
Perfluoro-n-heptanoic acid	PFHpA	2000
Perfluoro-n-octanoic acid	PFOA	2000
Perfluoro-n-nonanoic acid	PFNA	2000
Perfluoro-n-decanoic acid	PFDA	2000
Perfluoro-n-undecanoic acid	PFUdA	2000
Perfluoro-n-dodecanoic acid	PFDoA	2000
Perfluoro-n-tridecanoic acid	PFTrDA	2000
Perfluoro-n-tetradecanoic acid	PFTeDA	2000
2,3,3,3-Tetrafluoro-2-(1,1,2,2,3,3-heptafluoropropoxy)propanoic acid	HFPO-DA	2000
N-methylperfluoro-1-octanesulfonamidoacetic acid	N-MeFOSAA: linear isomer N-MeFOSAA: Σ branched isomers	1520 480
N-ethylperfluoro-1-octanesulfonamidoacetic acid	N-EtFOSAA: linear isomer N-EtFOSAA: Σ branched isomers	1550 450
Potassium perfluoro-1-butanesulfonate	L-PFBS	2000
Potassium perfluorohexanenesulfonate	PFHxSK: linear isomer PFHxSK: Σ branched isomers	1620 378
Potassium perfluorooctanenesulfonate	PFOSK: linear isomer PFOSK: Σ branched isomers	1580 422
Sodium dodecafluoro-3H-4,8-dioxanonanoate	NaDONA	1000
Potassium 9-chlorohexadecafluoro-3-oxanonane-1-sulfonate	9CI-PF3ONS	1000
Potassium 11-chloroeicosafuoro-3-oxaundecane-1-sulfonate	11CI-PF3OUdS	1000

用 EPA-537PDSL-R1	EPA Method 537 Analyte Primary Dilution Standard (linear isomers only)	1.2 ml
Perfluoro-n-hexanoic acid	PFHxA	2000
Perfluoro-n-heptanoic acid	PFHpA	2000
Perfluoro-n-octanoic acid	PFOA	2000
Perfluoro-n-nonanoic acid	PFNA	2000
Perfluoro-n-decanoic acid	PFDA	2000
Perfluoro-n-undecanoic acid	PFUdA	2000
Perfluoro-n-dodecanoic acid	PFDoA	2000
Perfluoro-n-tridecanoic acid	PFTrDA	2000
Perfluoro-n-tetradecanoic acid	PFTeDA	2000
2,3,3,3-Tetrafluoro-2-(1,1,2,2,3,3-heptafluoropropoxy)propanoic acid	HFPO-DA	2000
N-methylperfluoro-1-octanesulfonamidoacetic acid	N-MeFOSAA	2000
N-ethylperfluoro-1-octanesulfonamidoacetic acid	N-EtFOSAA	2000
Potassium perfluoro-1-butanesulfonate	L-PFBS	2000
Sodium perfluoro-1-hexanesulfonate	L-PFHxS	2000
Sodium perfluoro-1-octanesulfonate	L-PFOS	2000
Sodium dodecafluoro-3H-4,8-dioxanonanoate	NaDONA	2000
Potassium 9-chlorohexadecafluoro-3-oxanonane-1-sulfonate	9CI-PF3ONS	2000
Potassium 11-chloroeicosafuoro-3-oxaundecane-1-sulfonate	11CI-PF3OUdS	2000

NOTE: Listed concentrations are reported as the salt.

用 第一種特定化學物質

MIXED NATIVE PFAS: SOLUTION/MIXTURES

Catalogue Number	Product (methanol <1% water)	Qty/Conc
用 EPA-533PAR	EPA Method 533 Native Analyte Primary Dilution Standard (br/linear mix)	1.2 ml
Perfluoro-n-butanoic acid	PFBA	500
Perfluoro-n-pentanoic acid	PPPeA	500
Perfluoro-n-hexanoic acid	PFHxA	500
Perfluoro-n-heptanoic acid	PFHpA	500
Perfluoro-n-octanoic acid	PFOA	500
Perfluoro-n-nonanoic acid	PFNA	500
Perfluoro-n-decanoic acid	PFDA	500
Perfluoro-n-undecanoic acid	PFUdA	500
Perfluoro-n-dodecanoic acid	PFDoA	500
Perfluoro-n-tridecanoic acid	PFTrDA	500
Perfluoro-n-tetradecanoic acid	PFTeDA	500
2,3,3,3-Tetrafluoro-2-(1,1,2,2,3,3,3-heptafluoropropoxy)propanoic acid	HFPO-DA	500
Perfluoro-4-oxapentanoic acid	PF4OPeA	500
Perfluoro-5-oxahexanoic acid	PF5OHxA	500
Perfluoro-3,6-dioxaheptanoic acid	3,6-OPFHpA	500
Potassium perfluoro-1-butanesulfonate	L-PFBS	500
Sodium perfluoro-1-pentanesulfonate	L-PFPes	500
Potassium perfluorohexanenesulfonate	PFHxSK: linear isomer	406
	PFHxSK: Σ branched isomers	94.5
Sodium perfluoro-1-heptanesulfonate	L-PFHps	500
Potassium perfluoroctanenesulfonate	PFOSK: linear isomer	394
	PFOSK: Σ branched isomers	106
Sodium 1H,1H,2H,2H-perfluorohexane sulfonate	4:2FTS	500
Sodium 1H,1H,2H,2H-perfluorooctane sulfonate	6:2FTS	500
Sodium 1H,1H,2H,2H-perfluorodecane sulfonate	8:2FTS	500
Sodium dodecafluoro-3H-4,8-dioxanonanoate	NaDONA	500
Potassium 9-chlorohexadecafluoro-3-oxanonane-1-sulfonate	9CI-PF3ONS	500
Potassium 11-chloroeicosfluoro-3-oxaundecane-1-sulfonate	11CI-PF3OUdS	500
Potassium perfluoro(2-ethoxyethane)sulfonate	PFEESA	500
PFAC-MXF	Native Replacement PFAS Solution/Mixture	1.2 ml
2,3,3,3-Tetrafluoro-2-(1,1,2,2,3,3,3-heptafluoropropoxy)propanoic acid	HFPO-DA	2000
Sodium dodecafluoro-3H-4,8-dioxanonanoate	NaDONA	2000
Potassium 9-chlorohexadecafluoro-3-oxanonane-1-sulfonate	9CI-PF3ONS	2000
Potassium 11-chloroeicosfluoro-3-oxaundecane-1-sulfonate	11CI-PF3OUdS	2000
PFAC-MXG	Native Perfluoroether-Carboxylic Acids and -Sulfonate Solution/Mixture	1.2 ml
Perfluoro-4-oxapentanoic acid	PF4OPeA	2000
Perfluoro-5-oxahexanoic acid	PF5OHxA	2000
Perfluoro-3,6-dioxaheptanoic acid	3,6-OPFHpA	2000
Potassium perfluoro(2-ethoxyethane)sulfonate	PFEESA	2000

NOTE: Listed concentrations are reported as the salt.

用 第一種特定化学物質

MIXED NATIVE PFAS: SOLUTION/MIXTURES

Catalogue Number	Product (methanol <1% water)	Qty/Conc
NEW 用 ISO 21675-NSS	ISO 21675:2019 Native Stock Solution	1.2 ml
Perfluoro-n-butanoic acid	PFBA	100ng/ml
Perfluoro-n-pentanoic acid	PPPeA	100ng/ml
Perfluoro-n-hexanoic acid	PFHxA	100ng/ml
Perfluoro-n-heptanoic acid	PFHpA	100ng/ml
Perfluoro-n-octanoic acid	PFOA	100ng/ml
Perfluoro-n-nonanoic acid	PFNA	100ng/ml
Perfluoro-n-decanoic acid	PFDA	100ng/ml
Perfluoro-n-undecanoic acid	PFUdA	100ng/ml
Perfluoro-n-dodecanoic acid	PFDoA	100ng/ml
Perfluoro-n-tridecanoic acid	PFTrDA	100ng/ml
Perfluoro-n-tetradecanoic acid	PFTeDA	100ng/ml
Perfluoro-n-hexadecanoic acid	PFHxDA	100ng/ml
Perfluoro-n-octadecanoic acid	PFODA	100ng/ml
Perfluoro-1-octanesulfonamide	FOSA	100ng/ml
N-methylperfluoro-1-octanesulfonamide	N-MeFOSA	100ng/ml
N-ethylperfluoro-1-octanesulfonamide	N-EtFOSA	100ng/ml
N-methylperfluoro-1-octanesulfonamidoacetic acid	N-MeFOSAA	100ng/ml
N-ethylperfluoro-1-octanesulfonamidoacetic acid	N-EtFOSAA	100ng/ml
2H-Perfluoro-2-decenoic acid	FOUEA	100ng/ml
2,3,3,3-Tetrafluoro-2-(1,1,2,2,3,3,3-heptafluoropropoxy)-propanic acid	HFPO-DA	100ng/ml
Potassium perfluoro-1-butanesulfonate	L-PFBS	100ng/ml
Sodium perfluoro-1-hexanesulfonate	L-PFHxS	100ng/ml
Sodium perfluoro-1-heptanesulfonate	L-PFHpS	100ng/ml
Sodium perfluoro-1-octanesulfonate	L-PFOS	100ng/ml
Sodium perfluoro-1-decanesulfonate	L-PFDS	100ng/ml
Sodium 1H,1H,2H,2H-perfluorooctane sulfonate	6:2FTS	100ng/ml
Sodium 1H,1H,2H,2H-perfluorodecane sulfonate	8:2FTS	100ng/ml
Sodium dodecafluoro-3H-4,8-dioxanonanoate	NaDONA	100ng/ml
Potassium 9-chlorohexadecafluoro-3-oxanonane-1-sulfonate	9Cl-PF3ONS	100ng/ml
Sodium bis(1H,1H,2H,2H-perfluorodecyl)phosphate	8:2diPAP	100ng/ml

NOTE: Listed concentrations are reported as the salt.

用 第一種特定化學物質

MIXED MASS-LABELLED PFAS: SOLUTION/MIXTURES

Catalogue Number	Product (methanol <1% water)	Qty/Conc
EPA-537SS-R1	EPA Method 537 Surrogate Primary Dilution Standard (SUR PDS)	1.2 ml
Perfluoro-n-[1,2- ¹³ C ₂]hexanoic acid	MPFHxA	1000
Perfluoro-n-[1,2- ¹³ C ₂]decanoic acid	MPFDA	1000
2,3,3,3-Tetrafluoro-2-(1,1,2,2,3,3,3-heptafluoropropoxy)- ¹³ C ₃ -propanoic acid	M3HFPO-DA	1000
N-ethyl-d5-perfluoro-1-octanesulfonamidoacetic acid	d5-N-EtFOSAA	4000
EPA-537IS	EPA Method 537 Internal Standard Primary Dilution Standard (IS PDS)	1.2 ml
Perfluoro-n-[1,2- ¹³ C ₂]octanoic acid	M2PFOA	1000
N-methyl-d3-perfluoro-1-octanesulfonamidoacetic acid	d3-N-MeFOSAA	4000
Sodium perfluoro-1-[1,2,3,4- ¹³ C ₄]octanesulfonate	MPFOS	3000
Catalogue Number	Product (methanol/2% isopropanol, / <1% water)	Qty/Conc
MPFAC-24ES	Mass-Labelled PFAS Extraction Standard Solution	1.2 ml
Perfluoro-n-[1,2,3,4- ¹³ C ₄]butanoic acid	MPFBA	1000
Perfluoro-n-[1,2,3,4,5- ¹³ C ₅]pentanoic acid	M5PFPeA	1000
Perfluoro-n-[1,2,3,4,6- ¹³ C ₅]hexanoic acid	M5PFHxA	1000
Perfluoro-n-[1,2,3,4- ¹³ C ₄]heptanoic acid	M4PFHpA	1000
Perfluoro-n-[¹³ C ₈]octanoic acid	M8PFOA	1000
Perfluoro-n-[¹³ C ₉]nonanoic acid	M9PFNA	1000
Perfluoro-n-[1,2,3,4,5,6- ¹³ C ₆]decanoic acid	M6PFDA	1000
Perfluoro-n-[1,2,3,4,5,6,7- ¹³ C ₇]undecanoic acid	M7PFUdA	1000
Perfluoro-n-[1,2- ¹³ C ₂]dodecanoic acid	MPFDoA	1000
Perfluoro-n-[1,2- ¹³ C ₂]tetradecanoic acid	M2PFTeDA	1000
Perfluoro-1-[¹³ C ₈]octanesulfonamide	M8FOSA	1000
N-methyl-d3-perfluoro-1-octanesulfonamidoacetic acid	d3-N-MeFOSAA	1000
N-ethyl-d5-perfluoro-1-octanesulfonamidoacetic acid	d5-N-EtFOSAA	1000
Sodium perfluoro-1-[2,3,4- ¹³ C ₃]butanesulfonate	M3PFBS	1000
Sodium perfluoro-1-[1,2,3- ¹³ C ₃]hexanesulfonate	M3PFHxS	1000
Sodium perfluoro-[¹³ C ₈]octanesulfonate (>97%)	M8PFOS	1000
Sodium 1H,1H,2H,2H-perfluoro-1-[1,2- ¹³ C ₂]-hexane sulfonate	M2-4:2FTS	1000
Sodium 1H,1H,2H,2H-perfluoro-1-[1,2- ¹³ C ₂]-decane sulfonate	M2-6:2FTS	1000
N-ethyl-d5-perfluoro-1-octanesulfonamidoacetic acid	M2-8:2FTS	1000
Catalogue Number	Product (methanol <1% water)	Qty/Conc
EPA-533ES	EPA Method 533 Isotope Dilution Standard PDS	1.2 ml
Perfluoro-n-[1,2,3,4- ¹³ C ₄]butanoic acid	MPFBA	500
Perfluoro-n-[1,2,3,4,5- ¹³ C ₅]pentanoic acid	M5PFPeA	500
Perfluoro-n-[1,2,3,4,6- ¹³ C ₅]hexanoic acid	M5PFHxA	500
Perfluoro-n-[1,2,3,4- ¹³ C ₄]heptanoic acid	M4PFHpA	500
Perfluoro-n-[¹³ C ₈]octanoic acid	M8PFOA	500
Perfluoro-n-[¹³ C ₉]nonanoic acid	M9PFNA	500
Perfluoro-n-[1,2,3,4,5,6- ¹³ C ₆]decanoic acid	M6PFDA	500
Perfluoro-n-[1,2,3,4,5,6,7- ¹³ C ₇]undecanoic acid	M7PFUdA	500
Perfluoro-n-[1,2- ¹³ C ₂]dodecanoic acid	MPFDoA	500
2,3,3,3-Tetrafluoro-2-(1,1,2,2,3,3,3-heptafluoropropoxy)- ¹³ C ₃ -propanoic acid	M3HFPO-DA	500
Sodium perfluoro-1-[2,3,4- ¹³ C ₃]butanesulfonate	M3PFBS	500
Sodium perfluoro-1-[1,2,3- ¹³ C ₃]hexanesulfonate	M3PFHxS	500
Sodium perfluoro-[¹³ C ₈]octanesulfonate (>97%)	M8PFOS	500
Sodium 1H,1H,2H,2H-perfluoro-1-[1,2- ¹³ C ₂]-hexane sulfonate	M2-4:2FTS	2000
Sodium 1H,1H,2H,2H-perfluoro-1-[1,2- ¹³ C ₂]-decane sulfonate	M2-6:2FTS	2000
N-ethyl-d5-perfluoro-1-octanesulfonamidoacetic acid	M2-8:2FTS	2000
EPA-533IS	EPA Method 533 Isotope Performance Standard PDS	1.2 ml
Perfluoro-n-[2,3,4- ¹³ C ₃]butanoic acid	M3PFBA	1000
Perfluoro-n-[1,2- ¹³ C ₂]octanoic acid	M2PFOA	1000
Sodium perfluoro-1-[1,2,3,4- ¹³ C ₄]octanesulfonate	MPFOS	3000

NOTE: Listed concentrations are reported as the salt.

第一種特定化学物質

MIXED MASS-LABELLED PFAS: SOLUTION/MIXTURES


Catalogue Number	Product (methanol <1% water)	Qty/Conc
NEW 用 ISO 21675-LSS		
Perfluoro-n-[1,2,3,4- ¹³ C ₄]butanoic acid	ISO 21675:2019 Labelled Stock Solution	1.2 ml
Perfluoro-n-[1,2,3,4,5- ¹³ C ₅]pentanoic acid	MPFBA	100ng/ml
Perfluoro-n-[1,2,3,4,6- ¹³ C ₅]hexanoic acid	M5PFPeA	100ng/ml
Perfluoro-n-[1,2,3,4- ¹³ C ₄]heptanoic acid	M5PFHxA	100ng/ml
Perfluoro-n-[¹³ C ₈]octanoic acid	M4PFHpA	100ng/ml
Perfluoro-n-[¹³ C ₉]nonanoic acid	M8PFOA	100ng/ml
Perfluoro-n-[1,2,3,4,5,6- ¹³ C ₆]decanoic acid	M9PFNA	100ng/ml
Perfluoro-n-[1,2,3,4,5,6,7- ¹³ C ₇]undecanoic acid	M6PFDA	100ng/ml
Perfluoro-n-[1,2- ¹³ C ₂]dodecanoic acid	M7PFUdA	100ng/ml
Perfluoro-n-[1,2- ¹³ C ₂]tetradecanoic acid	MPFDa	100ng/ml
Perfluoro-n-[1,2- ¹³ C ₂]hexadecanoic acid	M2PFTeDA	100ng/ml
Perfluoro-1-[¹³ C ₈]octanesulfonamide	M2PFHxDA	100ng/ml
N-methyl-d3-perfluoro-1-octanesulfonamide	M8FOSA	100ng/ml
N-ethyl-d5-perfluoro-1-octanesulfonamide	d-N-MeFOSA	100ng/ml
N-methyl-d3-perfluoro-1-octanesulfonamidoacetic acid	d-N-EtFOSA	100ng/ml
N-ethyl-d5-perfluoro-1-octanesulfonamidoacetic acid	d3-N-MeFOSAA	100ng/ml
2H-Perfluoro-[1,2- ¹³ C ₂]-2-decanoic acid	d5-N-EtFOSAA	100ng/ml
2,3,3,3-Tetrafluoro-2-(1,1,2,2,3,3,3-heptafluoropropoxy)- ¹³ C ₃ -propanoic acid	MFOUEA	100ng/ml
Sodium perfluoro-1-[2,3,4- ¹³ C ₃]butanesulfonate	M3HFPO-DA	100ng/ml
Sodium perfluoro-1-hexane[¹⁸ O ₂]sulfonate	M3PFBS	100ng/ml
Sodium perfluoro-[¹³ C ₈]octanesulfonate	MPFHxS	100ng/ml
Sodium 1H,1H,2H,2H-perfluoro-1-[1,2- ¹³ C ₂]octane sulfonate	M8PFOS	100ng/ml
Sodium 1H,1H,2H,2H-perfluoro-1-[1,2- ¹³ C ₂]decane sulfonate	M2-6:2FTS	100ng/ml
Sodium bis(1H,1H,2H,2H-[1,2- ¹³ C ₂]perfluorodecyl)phosphate	M2-8:2FTS	100ng/ml
	M4-8:2diPAP	100ng/ml

NOTE: Listed concentrations are reported as the salt.


用 第一種特定化學物質

PFOS/PFOA ISOMERS


Catalogue Number	Product (methanol solution)	Qty/Conc	
P1MHpS	Perfluoro-1-methylheptane sulfonate	200µl	1.00 µg/ml
P3MHpS	Perfluoro-3-methylheptane sulfonate Perfluoro-3-methylheptanoic acid	200µl	1.00 µg/ml 1.90 µg/ml
P4MHpS	Perfluoro-4-methylheptane sulfonate Perfluoro-4-methylheptanoic acid	200µl	1.00 µg/ml 2.20 µg/ml
P5MHpS	Perfluoro-5-methylheptane sulfonate Perfluoro-5-methylheptanoic acid	200µl	1.00 µg/ml 1.96 µg/ml
P6MHpS	Perfluoro-6-methylheptane sulfonate Perfluoro-6-methylheptanoic acid	200µl	1.00 µg/ml 3.10 µg/ml
P55DMHxS	Perfluoro-5,5-dimethylhexane sulfonate Perfluoro-5,5-dimethyhexanoic acid	200µl	1.00 µg/ml 1.95 µg/ml
P45DMHxS	Perfluoro-4,5-dimethylhexane sulfonate Perfluoro-4,5-dimethylhexanoic acid Perfluoro-3,5-dimethylhexane sulfonate Perfluoro-3,5-dimethylhexanoic acid	200µl	1.00 µg/ml 1.22 µg/ml 0.50 µg/ml 0.60 µg/ml




Linear




Perfluoro-1-methyl-




Perfluoro-2-methyl-




Perfluoro-3-methyl-




Perfluoro-4-methyl-




Perfluoro-5-methyl-




Perfluoro-6-methyl-




Perfluoro-5,5-dimethyl-



Perfluoro-4,4-dimethyl-



Perfluoro-4,5-methyl-



Perfluoro-3,5-dimethyl-

NOTE: R=CO₂⁻ and CF₃SO₃⁻ except for Isomer 2 where R'=SO₃⁻ only

NATIVE PERFLUOROOCTANESULFONAMIDES (FOSAs)

Catalogue Number	Product	Qty/Conc	
FBSA-I	Perfluoro-1-butanesulfonamide in isopropanol	1.2 ml	50 µg/ml
FHxSA-I	Perfluoro-1-hexanesulfonamide in isopropanol	1.2 ml	50 µg/ml
FOSA-I	Perfluoro-1-octanesulfonamide in isopropanol	1.2 ml	50 µg/ml
FDSA-I	Perfluoro-1-decanesulfonamide in isopropanol	1.2 ml	50 µg/ml
N-MeFBSA-M	N-methylperfluoro-1-butanesulfonamide in methanol	1.2 ml	50 µg/ml
N-MeFOSA-M	N-methylperfluoro-1-octanesulfonamide in methanol	1.2 ml	50 µg/ml
N,N-Me2FOSA-M	N,N-dimethylperfluoro-1-octanesulfonamide in methanol	1.2 ml	50 µg/ml
N-EtFOSA-M	N-ethylperfluoro-1-octanesulfonamide in methanol	1.2 ml	50 µg/ml

MASS-LABELLED PERFLUOROOCTANESULFONAMIDES

Catalogue Number	Product
M8FOSA-I	 Perfluoro-1-[¹³ C ₈]octanesulfonamide 1.2 ml; 50 µg/ml ($\pm 2.5 \mu\text{g/ml}$); in isopropanol >99% linear; >99% ¹³ C ₈
d-N-MeFOSA-M	 N-methyl-d ₃ -perfluoro-1-octanesulfonamide 1.2 ml; 50 µg/ml ($\pm 2.5 \mu\text{g/ml}$); in methanol >99% linear; 98% ² H ₃
d-N-EtFOSA-M	 N-ethyl-d ₅ -perfluoro-1-octanesulfonamide 1.2 ml; 50 µg/ml ($\pm 2.5 \mu\text{g/ml}$); in methanol >99% linear; 98% ² H ₅

NATIVE PERFLUOROOCTANESULFONAMIDOETHANOLS (FOSEs)

Catalogue Number	Product (methanol solution)	Qty/Conc	
N-MeFOSE-M	2-(N-methylperfluoro-1-octanesulfonamido)-ethanol	1.2 ml	50 µg/ml
N-EtFOSE-M	2-(N-ethylperfluoro-1-octanesulfonamido)-ethanol	1.2 ml	50 µg/ml



MASS-LABELLED PERFLUOROOCTANESULFONAMIDOETHANOLS

Catalogue Number	Product
d7-N-MeFOSE-M	 2-(N-methyl-d ₃ -perfluoro-1-octane - sulfonamido)ethan-d ₄ -ol 1.2 ml; 50 µg/ml ($\pm 2.5 \mu\text{g/ml}$); in methanol >99% linear; 98% ² H ₇
d9-N-EtFOSE-M	 2-(N-ethyl-d ₅ -perfluoro-1-octane - sulfonamido)ethan-d ₄ -ol 1.2 ml; 50 µg/ml ($\pm 2.5 \mu\text{g/ml}$); in methanol >99% linear; 98% ² H ₉

NATIVE PERFLUOROOCTANESULFONAMIDOACETIC ACIDS (FOSAAs)

Catalogue Number	Product (methanol solution)	Qty/Conc	
FOSAA	Perfluoro-1-octanesulfonamidoacetic acid	1.2 ml	50 µg/ml
N-MeFOSAA	N-methylperfluoro-1-octanesulfonamidoacetic acid	1.2 ml	50 µg/ml
br-NMeFOSAA	N-methylperfluoro-1-octanesulfonamidoacetic acid isomeric mixture	1.2 ml	50 µg/ml
N-EtFOSAA	N-ethylperfluoro-1-octanesulfonamidoacetic acid	1.2 ml	50 µg/ml
br-NEtFOSAA	N-ethylperfluoro-1-octanesulfonamidoacetic acid isomeric mixture	1.2 ml	50 µg/ml







MASS-LABELLED PERFLUOROOCTANESULFONAMIDOACETIC ACIDS

Catalogue Number	Product
d3-N-MeFOSAA	 N-methyl-d ₃ -perfluoro-1-octane - sulfonamidoacetic acid 1.2 ml; 50 µg/ml (± 2.5 µg/ml); in methanol >99% linear; 98% ² H ₃
d5-N-EtFOSAA	 N-ethyl-d ₅ -perfluoro-1-octane - sulfonamidoacetic acid 1.2 ml; 50 µg/ml (± 2.5 µg/ml); in methanol >99% linear; 98% ² H ₅

NATIVE TELOMER ALCOHOLS (FTOHs)

Catalogue Number	Product (methanol solution)	Qty/Conc
FBET	2-Perfluorobutyl ethanol (4:2)	1.2 ml 50 µg/ml
5:2sFTOH	1-Perfluoropentyl ethanol (5:2 secondary)	1.2 ml 50 µg/ml
FHET	2-Perfluorohexyl ethanol (6:2)	1.2 ml 50 µg/ml
7:2sFTOH	1-Perfluoroheptyl ethanol (7:2 secondary)	1.2 ml 50 µg/ml
FOET	2-Perfluorooctyl ethanol (8:2)	1.2 ml 50 µg/ml
FDET	2-Perfluorodecyl ethanol (10:2)	1.2 ml 50 µg/ml

MASS-LABELLED TELOMER ALCOHOLS

Catalogue Number	Product
MFBET	 2-Perfluorobutyl-[1,1,2,2- ² H ₄]-ethanol (4:2) 1.2 ml; 50 µg/ml (± 2.5 µg/ml); in methanol ; > 98% 1,1,2,2- ² H ₄
MFHET	 2-Perfluorohexyl-[1,1- ² H ₂]-[1,2- ¹³ C ₂]-ethanol (6:2) 1.2 ml; 50 µg/ml (± 2.5 µg/ml); in methanol ; 98% 1,1- ² H ₂ >99% 1,2- ¹³ C ₂
M2FHET	 2-Perfluorohexyl- [1,2- ¹³ C ₂]-ethanol (6:2) 1.2 ml; 50 µg/ml (± 2.5 µg/ml); in methanol ; >99% 1,2- ¹³ C ₂
MFOET	 2-Perfluoroctyl-[1,1- ² H ₂]-[1,2- ¹³ C ₂]-ethanol (8:2) 1.2 ml; 50 µg/ml (± 2.5 µg/ml); in methanol ; 98% 1,1- ² H ₂ >99% 1,2- ¹³ C ₂
M2FOET	 2-Perfluoroctyl- [1,2- ¹³ C ₂]-ethanol (8:2) 1.2 ml; 50 µg/ml (± 2.5 µg/ml); in methanol ; >99% 1,2- ¹³ C ₂
MFDET	 2-Perfluorodecyl-[1,1- ² H ₂]-[1,2- ¹³ C ₂]-ethanol (10:2) 1.2 ml; 50 µg/ml (± 2.5 µg/ml); in methanol ; 98% 1,1- ² H ₂ >99% 1,2- ¹³ C ₂



NATIVE TELOMER ACIDS (FTAs)

Catalogue Number	Product (isopropanol solution)	Qty/Conc	
FHEA	2-Perfluorohexyl ethanoic acid (6:2)	1.2 ml	50 µg/ml
FOEA	2-Perfluoroctyl ethanoic acid (8:2)	1.2 ml	50 µg/ml
FDEA	2-Perfluorodecyl ethanoic acid (10:2)	1.2 ml	50 µg/ml
Catalogue Number	Product (methanol solution)		Qty/Conc
FPrPA	3-Perfluoropropyl propanoic acid (3:3)	1.2 ml	50 µg/ml
FPePA	3-Perfluoropentyl propanoic acid (5:3)	1.2 ml	50 µg/ml
FHpPA	3-Perfluoroheptyl propanoic acid (7:3)	1.2 ml	50 µg/ml

NATIVE TELOMER ACIDS: SOLUTION/MIXTURE

Catalogue Number	Product (isopropanol solution)	Qty/Conc
FTA-MXA	Native FTA Solution/Mixture	1.2 ml
2-Perfluorohexyl ethanoic acid (6:2)		2 µg/ml
2-Perfluoroctyl ethanoic acid (8:2)		2 µg/ml
2-Perfluorodecyl ethanoic acid (10:2)		2 µg/ml

MASS-LABELLED TELOMER ACIDS

Catalogue Number	Product
MFHEA	 2-Perfluorohexyl-[1,2- ¹³ C ₂]-ethanoic acid (6:2) 1.2 ml; 50 µg/ml (± 2.5 µg/ml); in isopropanol ; >99% ¹³ C ₂
MFOEA	 2-Perfluoroctyl-[1,2- ¹³ C ₂]-ethanoic acid (8:2) 1.2 ml; 50 µg/ml (± 2.5 µg/ml); in isopropanol ; >99% ¹³ C ₂
MFDEA	 2-Perfluorodecyl-[1,2- ¹³ C ₂]-ethanoic acid (10:2) 1.2 ml; 50 µg/ml (± 2.5 µg/ml); in isopropanol ; >99% ¹³ C ₂

MASS-LABELLED TELOMER ACIDS: SOLUTION/MIXTURE

Catalogue Number	Product (isopropanol solution)	Qty/Conc
MFTA-MXA	Mass-Labelled FTA Solution/Mixture	1.2 ml
2-Perfluorohexyl-[1,2- ¹³ C ₂]-ethanoic acid (6:2)		2 µg/ml
2-Perfluoroctyl-[1,2- ¹³ C ₂]-ethanoic acid (8:2)		2 µg/ml
2-Perfluorodecyl-[1,2- ¹³ C ₂]-ethanoic acid (10:2)		2 µg/ml

NATIVE UNSATURATED TELOMER ACIDS (FTUAs)

Catalogue Number	Product (isopropanol solution)	Qty/Conc
FHUEA	2H-Perfluoro-2-octenoic acid (6:2)	1.2 ml
FOUEA	2H-Perfluoro-2-decanoic acid (8:2)	1.2 ml
FDUEA	2H-Perfluoro-2-dodecanoic acid (10:2)	1.2 ml

MASS-LABELLED UNSATURATED TELOMER ACIDS

Catalogue Number	Product
MFHUEA	<p>2H-Perfluoro-[1,2-¹³C₂]-2-octenoic acid (6:2) 1.2 ml; 50 µg/ml (\pm2.5 µg/ml); in isopropanol ; > 99% ¹³C₂</p>
MFOUEA	<p>2H-Perfluoro-[1,2-¹³C₂]-2-decenoic acid (8:2) 1.2 ml; 50 µg/ml (\pm2.5 µg/ml); in isopropanol ; >99% ¹³C₂</p>
MFDUEA	<p>2H-Perfluoro-[1,2-¹³C₂]-2-dodecenoic acid (10:2) 1.2 ml; 50 µg/ml (\pm2.5 µg/ml); in isopropanol ; >99% ¹³C₂</p>

NATIVE CHLOROPERFLUOROETHER SULFONATES (CI-PFESAs)

Catalogue Number	Product (methanol solution)	Qty/Conc	
9CI-PF3ONS	Potassium 9-chlorohexadecafluoro-3-oxanonane-1-sulfonate	1.2 ml	50 µg/ml
11CI-PF3OUdS	Potassium 11-chloroeicosafafluoro-3-oxaundecane-1-sulfonate	1.2 ml	50 µg/ml

NATIVE POLYFLUOROETHER CARBOXYLATES (PFECAs)

Catalogue Number	Product (methanol solution)	Qty/Conc	
NaDONA	Sodium dodecafluoro-3H-4,8-dioxanonanoate	1.2 ml	50 µg/ml
PF4OPeA	Perfluoro-4-oxapentanoic acid (PFMPA)	1.2 ml	50 µg/ml
PF5OHxA	Perfluoro-5-oxahexanoic acid (PFMBA)	1.2 ml	50 µg/ml
3,6-OPFHpA	Perfluoro-3,6-dioxaheptanoic acid (NFDHA)	1.2 ml	50 µg/ml
P5MeODIOXOAc	Perfluoro([5-methoxy-1,3-dioxolan-4-yl]oxy)acetic acid (C6O4)	1.2 ml	50 µg/ml

NATIVE POLYFLUOROETHER SULFONATES (PFESAs)

Catalogue Number	Product (methanol solution)	Qty/Conc	
PFEESA	Potassium perfluoro(2-ethoxyethane)sulfonate	1.2 ml	50 µg/ml

NATIVE HEXAFLUOROPROPYLENE OXIDE DIMER ACID

Catalogue Number	Product
HFPO-DA	<p>2,3,3,3-Tetrafluoro-2-(1,1,2,2,3,3,3-heptafluoropropoxy)propanoic acid 1.2 ml; 50 µg/ml (\pm2.5 µg/ml); in methanol</p>

MASS-LABELLED HEXAFLUOROPROPYLENE OXIDE DIMER ACID

Catalogue Number	Product
M3HFPO-DA	<p>2,3,3,3-Tetrafluoro-2-(1,1,2,2,3,3,3-heptafluoropropoxy)-¹³C₃-propanoic acid 1.2 ml; 50 µg/ml (\pm2.5 µg/ml); in methanol</p>

NATIVE TELOMER SULFONATES(FTSs)

Catalogue Number	Product
4:2FTS	<p>Sodium 1H,1H,2H,2H,-perfluorohexane sulfonate (4:2) 1.2 ml x 50 µg/ml ($\pm 2.5\mu\text{g}/\text{ml}$); in methanol</p>
6:2FTS	<p>Sodium 1H,1H,2H,2H,-perfluoroctane sulfonate (6:2) 1.2 ml x 50 µg/ml ($\pm 2.5\mu\text{g}/\text{ml}$); in methanol</p>
8:2FTS	<p>Sodium 1H,1H,2H,2H,-perfluorodecane sulfonate (8:2) 1.2 ml x 50 µg/ml ($\pm 2.5\mu\text{g}/\text{ml}$); in methanol</p>
10:2FTS	<p>Sodium 1H,1H,2H,2H-perfluorododecane sulfonate (10:2) 1.2 ml x 50 µg/ml ($\pm 2.5\mu\text{g}/\text{ml}$); in methanol</p>

NOTE: Listed concentrations are reported as the salt.

MASS-LABELLED TELOMER SULFONATES (FTSs)


Catalogue Number	Product
M2-4:2FTS	<p>Sodium 1H,1H,2H,2H-perfluoro-1-[1,2- $^{13}\text{C}_2$]-hexane sulfonate (4:2) 1.2 ml x 50 µg/ml ($\pm 2.5\mu\text{g}/\text{ml}$); in methanol</p>
M2-6:2FTS	<p>Sodium 1H,1H,2H,2H-perfluoro-1-[1,2- $^{13}\text{C}_2$]-octane sulfonate (6:2) 1.2 ml x 50 µg/ml ($\pm 2.5\mu\text{g}/\text{ml}$); in methanol</p>
M2-8:2FTS	<p>Sodium 1H,1H,2H,2H-perfluoro-1-[1,2- $^{13}\text{C}_2$]-decane sulfonate (8:2) 1.2 ml x 50 µg/ml ($\pm 2.5\mu\text{g}/\text{ml}$); in methanol</p>

NOTE: Listed concentrations are reported as the salt.

NATIVE PERFLUOROALKYLPHOSPHONIC ACIDS (PFAPAs)

Catalogue Number	Product
PFHxPA	<p>Perfluorohexylphosphonic acid 1.2 ml; 50 µg/ml ($\pm 2.5\mu\text{g}/\text{ml}$); in methanol ;</p>
PFOPA	<p>Perfluoroctylphosphonic acid 1.2 ml; 50 µg/ml ($\pm 2.5\mu\text{g}/\text{ml}$); in methanol ;</p>
PFDPA	<p>Perfluorodecylphosphonic acid 1.2 ml; 50 µg/ml ($\pm 2.5\mu\text{g}/\text{ml}$); in methanol ;</p>
CI-PFHxPA	<p>6-Chloroperfluorohexylphosphonic acid 1.2 ml; 50 µg/ml ($\pm 2.5\mu\text{g}/\text{ml}$); in methanol ;</p>
CI-PFOPA	<p>8-Chloroperfluoroctylphosphonic acid 1.2 ml; 50 µg/ml ($\pm 2.5\mu\text{g}/\text{ml}$); in methanol ;</p>




NATIVE SODIUM PERFLUOROALKYL PHOSPHINATES (X:XPFPi)

Catalogue Number	Product
6:6 PFPI	 Sodium bis(perfluorohexyl)phosphinate 1.2 ml x 50 µg/ml (\pm 2.5µg/ml); in methanol
6:8 PFPI	 Sodium perfluorohexylfluorooctylphosphinate 1.2 ml x 50 µg/ml (\pm 2.5µg/ml); in methanol
8:8 PFPI	 Sodium bis(perfluorooctyl)phosphinate 1.2 ml x 50 µg/ml (\pm 2.5µg/ml); in methanol

NATIVE POLYFLUORINATED PHOSPHATE ESTERS (PAPs and diPAPs)

Catalogue Number	Product (methanol solution)	Qty/Conc	
6:2PAP	Sodium 1H,1H,2H,2H-perfluoroctylphosphate	1.2 ml	50 µg/ml
8:2PAP	Sodium 1H,1H,2H,2H-perfluorodecylphosphate	1.2 ml	50 µg/ml
6:2diPAP	Sodium bis(1H,1H,2H,2H-perfluoroctyl)phosphate	1.2 ml	50 µg/ml
6:2/8:2diPAP	Sodium (1H,1H,2H,2H-perfluoroctyl-1H,1H,2H,2H-perfluorodecyl)phosphate	1.2 ml	50 µg/ml
8:2diPAP	Sodium bis(1H,1H,2H,2H-perfluorodecyl)phosphate	1.2 ml	50 µg/ml

MASS-LABELLED POLYFLUORINATED PHOSPHATE ESTERS (M2-PAPs and M4-diPAPs)

Catalogue Number	Product
M2-6:2PAP	 Sodium 1H,1H,2H,2H-[1,2- ¹³ C ₂]perfluorooctylphosphate 1.2 ml; 50 µg/ml (\pm 2.5 µg/ml); in methanol ; >99% ¹³ C ₂
M2-8:2PAP	 Sodium 1H,1H,2H,2H-[1,2- ¹³ C ₂]perfluorodecylphosphate 1.2 ml; 50 µg/ml (\pm 2.5 µg/ml); in methanol ; >99% ¹³ C ₂
M4-6:2diPAP	 Sodium bis(1H,1H,2H,2H-[1,2- ¹³ C ₂]perfluorooctyl)-phosphate 1.2 ml; 50 µg/ml (\pm 2.5 µg/ml); in methanol ; >99% ¹³ C ₄
M4-8:2diPAP	 Sodium bis(1H,1H,2H,2H-[1,2- ¹³ C ₂]perfluorodecyl)-phosphate 1.2 ml; 50 µg/ml (\pm 2.5 µg/ml); in methanol ; >99% ¹³ C ₄

NATIVE POLYFLUORINATED PHOSPHATE ESTERS (SAmPAPs)

Catalogue Number	Product
SAmPAP	 Sodium 2-(N-ethylperfluorooctane-1-sulfonamido)ethyl phosphate 1.2 ml x 50 µg/ml ($\pm 2.5\mu\text{g}/\text{ml}$); in methanol
diSAmPAP	 Sodium bis-[2-(N-ethylperfluorooctane-1-sulfonamido)ethyl] phosphate 1.2 ml x 50 µg/ml ($\pm 2.5\mu\text{g}/\text{ml}$); in methanol

NATIVE TELOMER ACRYLATES (X:2FTAcr)

Catalogue Number	Product (isooctane solution)	Qty/Conc
8:2FTAcr	1H,1H,2H,2H-Perfluorodecyl acrylate	1.2 ml 50 µg/ml
10:2FTAcr	1H,1H,2H,2H-Perfluorododecyl acrylate (97%)	1.2 ml 50 µg/ml

NATIVE TELOMER ACETATES (X:2FTOAc)

Catalogue Number	Product (isooctane solution)	Qty/Conc
8:2FTOAc	1H,1H,2H,2H-Perfluorodecyl acetate	1.2 ml 50 µg/ml
10:2FTOAc	1H,1H,2H,2H-Perfluorododecyl acetate	1.2 ml 50 µg/ml

NATIVE AFFF PFAS

Catalogue Number	Product (methanol solution)	Qty/Conc
N-AP-FHxSA	1H,1H,2H,2H-Perfluorododecyl acetate	1.2 ml 50 µg/ml
N-TAmP-FHxSA	N-[3-(perfluoro-1-hexanesulfonamido)propan-1-yl]-N,N,N-trimethylammonium	1.2 ml 50 µg/ml
N-CMAmP-6:2FOSA	N-(carboxymethyl)-N,N-dimethyl-N-[3-(1H,1H,2H,2H-perfluoro-1-octanesulfonamido)propan-1-yl]ammonium (6:2FTAB)	1.2 ml 50 µg/ml
5:3FTB	2-[(4,4,5,5,6,6,7,7,8,8,8-Undecafluoroctyl)dimethylammonio]acetate	1.2 ml 50 µg/ml
5:1:2FTB	2-[(3,4,4,5,5,6,6,7,7,8,8,8-Dodecafluoroctyl)dimethylammonio]acetate	1.2 ml 50 µg/ml

用途確認書について

「第一種特定化学物質」に指定されている化合物については、製造、輸入、使用は原則として禁止されています。但し、試験研究用に使用する場合に限り、使用することが出来ます。この為、使用目的等を明記した「用途確認書」を弊社宛に提出頂き次第、販売させて頂きますのでご了承ください。

「化学物質の審査及び製造等の規制に関する法律(略称:化審法)」に指定されている第一種特定化学物質に該当する本カタログ掲載製品は、本文中に以下の略号で表示しています。

用 第一種特定化学物質

本カタログ掲載品中の第一種特定化学物質該当製品

• PFC-CVS-C	P1	• EU-5813-NSS	P7
• MPFAC-C-ES	P2	• MPFAC-MXA	P8
• MPFAC-C-IS	P2	• PFAC-24PAR	P8
• PFAC-MXC	P2	• PFAC-30PAR	P9
• L-PFOS	P3	• EPA-537PDS-R1	P10
• L-PFOSK	P3	• EPA-537PDSL-R1	P10
• PFS-MXA	P3	• EPA-533PAR	P11
• br-PFOSK	P3	• ISO 21675-NSS	P12
• T-PFOS	P3	• EPA-537IS	P13
• MPFOS	P4	• MPFAC-24ES	P13
• M8PFOS	P4	• EPA-533ES	P13
• PFAC-MXA	P7	• EPA-533IS	P13
• PFAC-MXB	P7	• ISO 21675-LSS	P14

ご購入をお考えの「第一種特定化学物質」該当製品がございましたら、下記までご連絡ください。
該当製品の「用途確認書ひな型」をお送り致します。



Wellington Laboratories Japan Inc.

株式会社 ウエリントン ラボラトリーズ ジャパン

〒177-0034 東京都練馬区富士見台1-22-12

TEL:03-5934-4184 E-mail:info@well-labs.co.jp

www.well-labs.co.jp



www.well-labs.co.jp

Wellington Laboratories Inc.

345 Southgate Drive, Guelph, ON, Canada N1G 3M5

www.well-labs.com