WELLINGTON Reporter

PRODUCT UPDATES FROM WELLINGTON LABORATORIES

November 2026

NEW PRODUCT

UK Drinking Water Inspectorate (DWI) Native PFAS Solution/Mixture

UK-DWI-MXA

Variation in analyte target lists for per- and polyfluoroalkyl substances (PFAS) between different geographical areas is a common issue and reflects a rapidly evolving regulatory environment. Diverse regional guidelines necessitate the tailoring of PFAS standards in order to meet local requirements. In the United Kingdom, the Drinking Water Inspectorate (DWI) regulates PFAS concentrations in public drinking water. In order to ensure that the PFAS levels do not pose a threat to human health, they have expanded their list of target compounds to include analytes of local concern. In response to customer inquiries, **Wellington** has developed a certified reference standard solution/mixture (**UK-DWI-MXA**) that can be used in conjunction with **EPA-1633MXA** to meet the requirements of this regional directive.

Wellington is committed to working with our customers to support PFAS analyses globally. If you have a product suggestion that can improve the applicability of our solution/mixtures in your area, please submit it for consideration via our webform (https://well-labs.com/newproducts/product-suggestion/) or contact us directly at info@well-labs.com.

	Product Code	Description (methanol solution)		Units/Conc
NEW	**UK-DWI-MXA**	UK DWI Native PFAS Solution/Mixture		1.2 mL
	Perfluoro-n-hexadecanoic acid		PFHxDA	2000 ng/mL
	Perfluoro-n-octadecanoic acid		PFODA	2000 ng/mL
	Hexafluoropropylene oxide trimer acid		HFPO-TrA	2000 ng/mL
	Perfluoro-1-butanesulfonamide		FBSA	2000 ng/mL
	Perfluoro-1-hexanesulfonamide		FHxSA	2000 ng/mL
	N-[3-(Dimethylamino)propan-1-yl]perfluoro-1-hexanesulfonamide		N-AP-FHxSA	2000 ng/mL
	6:2 Fluorotelomer sulfonamide alkylbetaine (6:2 FTAB)		N-CMAmp-6:2FOSA	2000 ng/mL
	Sodium perfluoro-1-undecanesulfonate		L-PFUdS	2000 ng/mL
	Sodium perfluoro-1-tridecanesulfonate		L-PFTrDS	2000 ng/mL
	Potassium perfluoro-4-ethylcyclohexanesulfonate (*cis*/*trans* isomer mixture)		PFECHS	2000 ng/mL
	Ammonium perfluoro-3,6-dioxaoctanoate		3,6-OPFOA	2000 ng/mL
	Sodium 1H,1H,2H,2H-perfluorododecanesulfonate		10:2FTS	2000 ng/mL
	Sodium bis(1H,1H,2H,2H-perfluorooctyl) phosphate		6:2diPAP	2000 ng/mL
	Sodium bis(1H,1H,2H,2H-perfluorodecyl) phosphate		8:2diPAP	2000 ng/mL

NOTE: Concentrations are listed for the specified salt when a counterion is indicated.

Please contact your local distributor or info@well-labs.com for pricing and delivery.

Visit our website (www.well-labs.com) for a complete listing of new products.

Product Code	Description (methanol solution)	Units
EPA-1633MXA	U.S. EPA Method 1633 Native PFAS Stock Solution/Mixture	1.2 mL

		EPA-1633MXA (ng/mL)
Perfluoro-n-butanoic acid	PFBA	1000
Perfluoro-n-pentanoic acid	PFPeA	1000
Perfluoro-n-hexanoic acid	PFHxA	1000
Perfluoro-n-heptanoic acid	PFHpA	1000
Perfluorooctanoic acid (linear and branched isomers)	br-PFOA	1000
Perfluorononanoic acid (linear and branched isomers)	br-PFNA	1000
Perfluoro-n-decanoic acid	PFDA	1000
Perfluoro-n-undecanoic acid	PFUdA	1000
Perfluoro-n-dodecanoic acid	PFDoA	1000
Perfluoro-n-tridecanoic acid	PFTrDA	1000
Perfluoro-n-tetradecanoic acid	PFTeDA	1000
Perfluorooctanesulfonamide (linear and branched isomers)	br-FOSA	1000
N-Methylperfluorooctanesulfonamide (linear and branched isomers)	br-NMeFOSA	1000
N-Ethylperfluorooctanesulfonamide (linear and branched isomers)	br-NEtFOSA	1000
N-Methylperfluorooctanesulfonamidoacetic acid (linear and branched isomers)	br-NMeFOSAA	1000
N-Ethylperfluorooctanesulfonamidoacetic acid (linear and branched isomers)	br-NEtFOSAA	1000
2-(N-Methylperfluorooctanesulfonamido)ethanol (linear and branched isomers)	br-NMeFOSE	1000
2-(N-Ethylperfluorooctanesulfonamido)ethanol (linear and branched isomers)	br-NEtFOSE	1000
rac-2,3,3,3-Tetrafluoro-2-(1,1,2,2,3,3,3-heptafluoropropoxy)propanoic acid	HFPO-DA	1000
Perfluoro-4-oxapentanoic acid	PF4OPeA	1000
Perfluoro-5-oxahexanoic acid	PF5OHxA	1000
Perfluoro-3,6-dioxaheptanoic acid	3,6-OPFHpA	1000
3-Perfluoropropyl propanoic acid	FPrPA	1000
3-Perfluoropentyl propanoic acid	FPePA	1000
3-Perfluoroheptyl propanoic acid	FHpPA	1000
Potassium perfluoro-1-butanesulfonate	L-PFBS	1000
Sodium perfluoro-1-pentanesulfonate	L-PFPeS	1000
Potassium perfluorohexanesulfonate (linear and branched isomers)	br-PFHxSK	1000
Sodium perfluoro-1-heptanesulfonate	L-PFHpS	1000
Potassium perfluorooctanesulfonate (linear and branched isomers)	br-PFOSK	1000
Sodium perfluoro-1-nonanesulfonate	L-PFNS	1000
Sodium perfluoro-1-decanesulfonate	L-PFDS	1000
Sodium perfluoro-1-dodecanesulfonate	L-PFDoS	1000
Sodium 1H,1H,2H,2H-perfluorohexanesulfonate	4:2FTS	1000
Sodium 1H,1H,2H,2H-perfluorooctanesulfonate	6:2FTS	1000
Sodium 1H,1H,2H,2H-perfluorodecanesulfonate	8:2FTS	1000
Sodium *rac*-dodecafluoro-3H-4,8-dioxanonanoate	NaDONA	1000
Potassium 9-chlorohexadecafluoro-3-oxanonane-1-sulfonate	9Cl-PF3ONS	1000
Potassium 11-chloroeicosafluoro-3-oxaundecane-1-sulfonate	11Cl-PF3OUdS	1000
Potassium perfluoro(2-ethoxyethane)sulfonate	PFEESA	1000

NOTE: Concentrations are listed for the specified salt when a counterion is indicated.